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English

1. Product description


The company's oil-free air compressor uses a high-quality, high-flow oil-free piston compressor as its power source. This provides a stable source of oil-free air, prevents oil particles from entering the human body, and avoids damage to end equipment caused by oil in the piping system. The components are selected for high quality. These products are widely used in the medical field, food processing, household applications, environmental protection, biopharmaceutical applications, and other sectors. The machine is fully designed in line with the requirements of dental applications: high air output, low noise, a dry and clean air source, smooth operation, and automatic control. During operation, the compressor starts or stops automatically when tank pressure reaches the preset minimum or maximum value. The compressor is especially suitable as a gas source for dental treatment and is an ideal choice for small and medium-sized dental hospitals and private clinics.

2. Product specifications

Model	Rated voltage (V)	Frequency (Hz)	Rated power (W)	Capacity (L/min)	Noise (dB)	Working pressure (MPa)	Max. outlet pressure (MPa)	Tank (L)	Dimensions (L×W×H) mm	Note
EWS06	220-240 / 110-220	50 / 60	480	89	48	0,8	0,85	6	410×180×470	
EWS24	220-240 / 110-220	50 / 60	680	116	50	0,8	0,85	24	530×300×520	
EWS30	220-240 / 110-220	50 / 60	800	150	55	0,8	0,85	30	600×350×570	
EWS40	220-240 / 110-220	50 / 60	1200	196	62	0,8	0,85	40	670×400×660	

3. Exploded view and parts list

Original exploded view from the manual.

1	cover	17	valve	33	stator	49	nut
2	left fan	18	valve	34	bearing6203-2Z	50	wheel
3	leftcrankcase	19	valve plate	35	Abjustingwasher	51	bolt
4	connectingrod	20	washer	36	vibrationcolumn	52	drain cock
5	piston cup	21	right fan	37	spring	53	bolt
6	Binder plate	22	shaft seal	38	air filter	54	washer
7	bolt	23	bearing6908-2Z	39	elbow	55	cushion foot
8	cylinder	24	bolt	40	Deflationvalve	56	washer
9	obturatorring	25	bolt	41	regulator	57	nut
10	obturatorring	26	right canckcase	42	40gauge	58	pipe
11	cylinder head	27	protectionring	43	50gauge	59	actuated valve
12	bolt	28	bolt	44	connect	60	checkvalve
13	obturatorring	29	bolt	45	pressureswitch	61	plug
14	connectingpipe	30	spring	46	safety valve	62	capacitor
15	bolt	31	rotator	47	handle	63	capacitor cover
16	limitedblock	32	nut	48	tank		

No.	Part	No.	Part
1	cover	2	left fan
3	left crankcase	4	connecting rod
5	piston cup	6	binder plate
7	bolt	8	cylinder
9	sealing ring	10	sealing ring
11	cylinder head	12	bolt
13	sealing ring	14	connecting pipe
15	bolt	16	limit block
17	valve	18	valve
19	valve plate	20	washer
21	right fan	22	shaft seal
23	bearing 6908-2Z	24	bolt
25	bolt	26	right crankcase
27	protection ring	28	bolt
29	bolt	30	spring
31	rotator	32	nut
33	stator	34	bearing 6203-2Z
35	adjusting washer	36	vibration column
37	spring	38	air filter
39	elbow	40	deflation valve
41	regulator	42	40 gauge
43	50 gauge	44	connector
45	pressure switch	46	safety valve
47	handle	48	tank
49	nut	50	wheel
51	bolt	52	drain cock
53	bolt	54	washer
55	cushion foot	56	washer
57	nut	58	pipe
59	actuated valve	60	check valve
61	plug	62	capacitor
63	capacitor cover		

4. Use and operating instructions

4.1 During use, place the air compressor on a level, dry, clean, and well-ventilated indoor surface. Pay special attention to keeping it away from dust. If these conditions cannot be met, additional filtration must be provided.

4.2 Before use, check whether the gas circuit of the compressor is blocked and whether any component shows abnormal conditions. Then plug the power cord into a suitable, properly grounded outlet. The motor starts, the solenoid valve operates, the vent line closes, and the unit begins to run.

4.3 When the pressure gauge indicates that the preset maximum pressure has been reached, the pressure control automatically disconnects and the compressor stops. If the compressor did not stop at maximum pressure because of a power interruption, it will resume operation after power is restored once the tank pressure drops to the lower pressure-control point.

4.4 Connect the valve outlet to the pipeline of the compressed-air equipment. When the outlet valves are open, compressed air can be supplied directly to the equipment. When the pressure gauge indicates that pressure has fallen to the preset minimum pressure, the pressure control automatically closes again and the compressor restarts. The minimum and maximum pressure values are factory-set for each model; users generally cannot adjust these settings themselves.

4.5 If the supply voltage is more than 10% below the rated voltage, the compressor must not be forced into use; otherwise the electrical system may be damaged. In that case, use a voltage stabilizer or booster so that the supply voltage remains within the required range.

4.6 During operation, a small amount of moisture condenses in the air tank of the compressor. It must be drained regularly, generally about once a week. Before draining, switch off the power, open the outlet valves to release the air from the tank, and then loosen the bottom drain plug or seal. Reinstall it after the water has been drained.

4.7 Because the intake air may contain dust, the silencer should be cleaned once a week. Before cleaning, switch off the power, remove the silencer, remove the end cap, take out the sponge inside, clean it thoroughly with clean water, let it dry, and reinstall it. The sponge may be replaced if necessary.

4.8 A pressure regulator with an adjusting screw is installed on the air line and can be used to regulate pressure. The manufacturer adjusts this component before delivery. Users must not change that setting at their own discretion, otherwise the design limits of the equipment may be exceeded, leading to damage or abnormal operation.

4.9 The compressor is equipped with thermal protection and overcurrent protection. When the motor overheats, the thermal protector disconnects automatically; when the motor stalls or the voltage is too high, the overcurrent protector disconnects automatically. Both protections reset automatically. During inspection, therefore, ensure that the circuit is unobstructed and can flow freely.

5. Maintenance and protection

5.1 Do not touch any component of the appliance while the compressor is operating.

5.2 This is an oil-free compressor. Adding oil or lubricant is prohibited.

5.3 The ambient temperature for operating this compressor is 5°C to 40°C.

5.4 Do not wet the complete appliance while cleaning the compressor. Use a soft dry cloth to remove dust and avoid damaging the components.

5.5 After use, switch off the power and shut off the air supply.

6. Fault analysis and remedies

No.	Fault characteristic	Possible cause	Remedy
1	The compressor does not start; the electric motor makes no humming sound.	There is no power at the plug.	Switch on the power supply.
		The overcurrent protector or pressure switch is not connected.	Check whether the connections of the overcurrent protector and pressure controller are properly connected.
2	The compressor stops before the preset maximum pressure is reached, and air can	The set value of the pressure switch is too low.	Adjust the pressure switch to the required value.


No.	Fault characteristic	Possible cause	Remedy
	be heard venting through the solenoid valve.		
		The overcurrent protector is faulty.	Replace the overcurrent protector.
3	The compressor does not start normally; the motor turns but stops again immediately and a humming sound is heard for 10 to 30 seconds.	The voltage is too low.	Increase or stabilize the supply voltage. Use the compressor again only after the voltage has returned to normal.
		There is still air in the buffer tank, or the plug does not make good contact, or the power interruption lasted less than five seconds.	Unplug the unit. Replace the plug if necessary to ensure proper contact. Check whether the safety valve is defective. Procedure: unplug the unit, release the air from the tank, plug it in again, and check whether the solenoid valve vents when the required pressure is reached.
		The capacitor is damaged or its capacity has decreased.	Replace the capacitor.
		The compressor itself is faulty.	Contact the supplier.
4	The compressor keeps running, but the pressure does not reach the cut-off pressure or takes longer than normal to reach it.	There is leakage in the compressor or piping system.	Check all joints for air leakage, for example with soapy water.
		The compressor pump is faulty.	Open the cylinder head and check whether the valve has carbon deposits. Check whether the valve plate is deformed and whether the leather cup is worn. Contact the supplier if necessary.
5	After the compressor stops and the vent plug is switched off, atmospheric pressure drops.	Air leakage at the tank connections.	Apply sealant and tighten the connections.
		The holding valve does not seal airtight; leakage or reverse flow is present.	Clean the valve plug thoroughly. Replace the holding valve or valve plug if necessary.
6	The compressor does not stop when the preset maximum pressure is exceeded.	The pressure controller is faulty.	Adjust or replace the pressure controller.
7	Within the compressor working-pressure range, the safety valve vents air.	The spring force of the safety valve has weakened or the adjusting nut is loose.	Readjust the safety valve to the correct working pressure and then lock the nut.

7. Gas circuit diagram

Original gas-circuit diagram from the manual.

5	plug, the atmospheric pressure went down.	Holding valve has no airproof, steam leakage reversing.	Clear, clean valve plug. Replace holding valve or valve plug.
6	The compressor can't stop when is over the set maximum pressure.	Pressure controller broken.	Adjust or replace the pressure controller.
7	In the scope of compressor's work pressure, the safety valve exhaust.	Safety valve spring force weakened, the adjustment nut loose	Adjust safety valve adjust the nut to the safety valve working pressure exhaust, and locking nut.

7. Gas circuit elementary diagram

Translation of text in image

Inlet filter muffler: Inlet filter muffler

air compressor mainframe: air compressor mainframe

holding valve: holding valve

pressure gauge: pressure gauge

Gas tank: gas tank / air tank

pressure switch: pressure switch

Disconnecting valve: disconnecting valve

safety valve: safety valve

drain valve: drain valve